

Sepsis MPage - Data Display Specification

Overview

This document specifies the data display logic for the Sepsis MPage application, with particular emphasis on intervention column logic and clinical decision support features.

Table of Contents

- 1. [Column Structure](#)
- 2. [Intervention Display Logic](#)
- 3. [Clinical Decision Logic](#)
- 4. [Visual Styling Rules](#)
- 5. [Data Integration](#)
- 6. [SEP-1 Bundle Compliance](#)

Column Structure

Complete Column Reference

Index	Column Name	Data Field	Width	Type	Description
0	Patient	PATIENT_NAME	120px	text	Patient name (ZZZTEST prefix for test data)
1	Unit	UNIT	50px	text	Hospital unit (ER, ICU, 3W, 2N)
2	Room/Bed	ROOM_BED	80px	text	Room and bed number
3	Alert	ALERT_TYPE	50px	html	Alert severity indicator with clickable icon
4	Priority	SEPSIS_PRIORITY	60px	text	Sepsis priority level (High/Medium)
5	Time Zero	LAST_SEPSIS_SCREEN	140px	text	Sepsis identification timestamp
6	Timer	SEPSIS_TIMER	100px	text	Elapsed time since Time Zero
7	PowerPlan	POWERPLAN_ORDERED	75px	text	Sepsis PowerPlan order status
8	① Lac	LACTATE_ORDERED	50px	text	Initial lactate order status
9	① Result	LACTATE_RESULT	60px	numeric	Initial lactate result value
10	Cultures	BLOOD_CULTURES_ORDERED	65px	text	Blood culture order status
11	Abx	ANTIBIOTICS_ORDERED	45px	text	Antibiotics order status
12	② Lac	REPEAT_LACTATE_ORDERED	50px	text	Repeat lactate order status
13	Fluids	SEPSIS_FLUID_ORDERED	60px	text	IV fluid order status
14	Vol Doc	FLUID_VOLUME_DOCUMENTED	60px	text	Fluid volume documentation status
15	Perfusion	PERFUSION_ASSESSMENT	70px	text	Perfusion assessment status
16	Pressors	PRESSORS_ORDERED	65px	text	Vasopressor order status
17	Timeline	SHOW_TIMELINE	60px	html	Timeline visualization (Currently Hidden)

Intervention Display Logic

Standard Intervention Columns

Applies to: PowerPlan, ① Lac, Cultures, Abx, ② Lac, Fluids, Vol Doc

Data Values and Display

Data Value	Display	Icon	Color	Font Weight	CSS Class
"Y"	✓	fas fa-check	Green (#16a34a)	Bold	text-green-600 font-bold
"N"	✗	fas fa-times	Red (#dc2626)	Bold	text-red-600 font-bold
"Pend"		fas fa-hourglass-half	Orange (#f59e0b)	Normal	text-yellow-600 italic
null, "", missing	--	Text	Gray (#9ca3af)	Normal	N/A

Implementation Code

```
if (value === "Y") {
  TD.innerHTML = '<i class="fas fa-check"></i>';
  TD.style.color = "#16a34a";
  TD.style.fontWeight = "bold";
} else if (value === "N") {
  TD.innerHTML = '<i class="fas fa-times"></i>';
  TD.style.color = "#dc2626";
  TD.style.fontWeight = "bold";
} else if (value === "Pend") {
  TD.innerHTML = '<i class="fas fa-hourglass-half"></i>';
  TD.style.color = "#f59e0b";
  TD.style.fontWeight = "normal";
} else {
  TD.innerHTML = "--";
  TD.style.color = "#9ca3af";
}
```

Lactate Result Column (① Result) - Special Logic

Display Rules

1. When lactate is ordered and has result:
`javascript if (LACTATE_ORDERED === "Y" && LACTATE_RESULT exists) { Display: LACTATE_RESULT.toFixed(1) // e.g., "4.5" }`
2. Critical value highlighting (≥ 4.0 mmol/L):
`javascript if (LACTATE_RESULT >= 4.0) { Background: Pink (#fee2e2) Text Color: Red (#dc2626) Font Weight: Bold }`
3. No result available: `javascript Display: "--" Color: Gray (#9ca3af)`

Clinical Significance

- **Normal:** < 2.0 mmol/L
- **Elevated:** 2.0 - 3.9 mmol/L (triggers repeat lactate requirement)
- **Critical:** ≥ 4.0 mmol/L (indicates severe tissue hypoperfusion, triggers perfusion/pressor protocols)

Perfusion & Pressors Columns - Conditional Logic

Special Display Rules

These columns have conditional display based on lactate results:

```
// Columns 15 (Perfusion) and 16 (Pressors)
if (LACTATE_RESULT < 4.0 || LACTATE_RESULT === null || LACTATE_RESULT === undefined) {
  Display: "N/A"
  Color: Gray (#9ca3af)
  Font Weight: Normal
  Rationale: Not clinically indicated
} else {
  // Display standard Y/N/Pend based on actual data
  Apply standard intervention display logic
}
```

Clinical Rationale

- Perfusion reassessment and vasopressors are only indicated when lactate ≥ 4.0 mmol/L
- "N/A" prevents inappropriate alerts for interventions that aren't clinically warranted
- Follows SEP-1 guidelines for shock management

Clinical Decision Logic

Alert Types and Criteria

1. SIRS (Systemic Inflammatory Response Syndrome)

- **Icon:** Yellow circle (fas fa-circle , #eab308)
- **Criteria:** 3 or more of the following:
 - Heart Rate > 90 bpm
 - Respiratory Rate > 20/min
 - Temperature > 38.3°C or < 36°C
 - WBC < 4 or > 12 x10⁹/L
 - Glucose 7.8-11.1 mmol/L (in non-diabetic)

2. Sepsis (SIRS + Organ Dysfunction)

- **Icon:** Orange circle (fas fa-circle , #f97316)
- **Criteria:** SIRS + one or more organ dysfunction markers:
 - Lactate > 2.0 mmol/L
 - Creatinine > 176.8 μmol/L
 - Bilirubin 35.9-171 μmol/L
 - Platelets < 100 x10⁹/L
 - INR > 1.5

3. Septic Shock (SIRS + Cardiovascular Dysfunction)

- **Icon:** Red circle (fas fa-circle , #dc2626)
- **Criteria:** SIRS + cardiovascular dysfunction:
 - SBP < 90 mmHg
 - MAP < 65 mmHg
 - Lactate ≥ 4.0 mmol/L

Timer Display and Bundle Compliance

Time Format

- **Under 1 hour:** Display as `Xm` (e.g., "45m")
- **1 hour or more:** Display as `Xh Ym` (e.g., "2h 15m")

Visual Indicators by Time Elapsed

Time Range	Status	Text Color	Background	Meaning
0-2 hours	On track	Default	None	Within optimal window
2-3 hours	Warning	Yellow-brown (#ca8a04)	None	Approaching 3-hour deadline
3-6 hours	Missed 3-hour	Orange (#ea580c)	Light orange (#fed7aa)	3-hour bundle missed
>6 hours	Critical	Red (#dc2626)	Light red (#fee2e2)	Both bundles missed

Visual Styling Rules

Color Palette

Purpose	Color	Hex Code	Usage
Success/Complete	Green	#16a34a	Completed interventions (Y)
Error/Not Done	Red	#dc2626	Not ordered interventions (N)
Warning/Pending	Orange	#f59e0b	Pending interventions
Inactive/N/A	Gray	#9ca3af	Not applicable or no data
Critical Background	Pink	#fee2e2	Lactate ≥ 4.0, Timer >6h
Warning Background	Light Orange	#fed7aa	Timer 3-6h

Font Awesome Icons

- **Version:** 6.5.1
- **Icons Used:**
 - `fas fa-check` - Completed/Yes
 - `fas fa-times` - Not ordered/No
 - `fas fa-hourglass-half` - Pending
 - `fas fa-circle` - Alert severity indicators

Data Integration

Expected JSON Structure

```

{
  "REC": {
    "CCL_SCRIPT": "1_MP_SEPSIS_HANDSONTABLE",
    "SEPSIS": [
      {
        "PATIENT_NAME": "ZZZTEST, Name",
        "UNIT": "ICU",
        "ROOM_BED": "201-A",
        "ALERT_TYPE": "SIRS",
        "ALERT_DETAILS": {
          "HR": 96,
          "RR": 22,
          "TEMP": 38.1,
          "WBC": 13.2,
          "LACTATE": 2.3
        },
        "SEPSIS_PRIORITY": "High",
        "LAST_SEPSIS_SCREEN": "08/13/25 16:06 [Yes]",
        "SEPSIS_TIMER": "00:24",
        "POWERPLAN_ORDERED": "Y",
        "LACTATE_ORDERED": "Y",
        "LACTATE_RESULT": 4.5,
        "ANTIBIOTICS_ORDERED": "Pend",
        "BLOOD_CULTURES_ORDERED": "Y",
        "SEPSIS_FLUID_ORDERED": "N",
        "REPEAT_LACTATE_ORDERED": "N",
        "PERFUSION_ASSESSMENT": "N",
        "PRESSORS_ORDERED": "N"
      }
    ]
  }
}

```

Data Validation Rules

1. **String fields:** Accept "Y", "N", "Pend", or null/empty
2. **Numeric fields:** LACTATE_RESULT accepts float values
3. **Missing fields:** Default to "---" display
4. **Invalid values:** Numeric values in Y/N/Pend fields display as "---"
5. **Timer format:** Expected as "HH:MM" string

SEP-1 Bundle Compliance

3-Hour Bundle (Must complete within 3 hours)

1. **Measure lactate level** (① Lac + ① Result)
2. **Obtain blood cultures** before antibiotics (Cultures)
3. **Administer broad-spectrum antibiotics** (Abx)
4. **Begin rapid administration of 30mL/kg crystalloid** for hypotension or lactate ≥ 4 (Fluids)

6-Hour Bundle (Must complete within 6 hours)

1. **Re-measure lactate** if initially elevated >2 (② Lac)
2. **Apply vasopressors** if hypotensive after fluid resuscitation (Pressors)
3. **Reassess volume status and tissue perfusion** (Perfusion)

Bundle Compliance Indicators

- **Green background zones:** Interventions completed within time windows
- **Timer colors:** Visual cues for bundle deadline status

- **Column ordering:** Follows SEP-1 intervention sequence

Implementation Notes

Performance Considerations

- Handsontable renders only visible cells (virtual rendering)
- AfterRenderer called on scroll/update - keep logic lightweight
- Font Awesome icons loaded via CDN for consistency

Accessibility

- Color coding supplemented with icons (not color alone)
- Text alternatives for all visual indicators
- Clickable elements have cursor pointer and tooltips

Future Enhancements

- Timeline visualization (currently hidden, column 17)
- Real-time timer updates
- CCL backend integration for live data
- Export functionality for compliance reporting

=====

Last Updated: August 20, 2025 Version: 1.0